

ABBREVIATED IVP PROTOCOL

This checklist can be used to ensure that each step in the IVP procedure is completed successfully.

DAY -2

Preparation of Media

- saline
- OMM
- OCM

Check incubators for accuracy of temperature, O₂ and CO₂ readings

DAY -1

Collection of ovaries from slaughterhouse

Preparation for oocyte collection

- Place OCM and saline in the oven
- Prepare OMM microdrops and cover them with mineral oil
- Place plates in the incubator
- Set up for oocyte collection
 - Scalpel
 - Scalpel blades
 - Gloves
 - 400 ml beaker
 - container to discard ovaries
 - bench paper to cover surface

Harvesting oocyte-cumulus complexes from ovaries

- Add 100 ml of OCM to beaker
- Slice ovaries
- Swirl ovaries in the beaker
- Pour OCM and oocyte mixture into 50 ml sterile centrifuge tubes and place into water bath

Collection of cumulus oocyte complexes (COCs) from centrifuge tubes (5 min. sedimentation time)

- Set up for oocyte collection
 - 100 µm cell strainer
 - sterile transfer pipets
 - 10 ml airtite syringe
 - 18 gauge needle
 - Grid plate
- Set up cell strainer over 50 ml beaker

- Using transfer pipet, suck up pellet and place in cell strainer
- Flip cell strainer over into a grid plate on a slide warmer
- Fill needle/syringe with OCM and rinse the debris from the strainer into the grid dish

Searching for cumulus oocyte complexes

- X-plate
- Dissecting microscope
- Searching instrument (microdispensor, wiretrol, etc.)
- Slide warmer
- Transfer cumulus oocyte complexes to X-plate and rinse two times

Oocyte maturation

- After the last rinsing place cleaned cumulus oocytes complexes (10/drop) into a 50 μ l microdrop of pre-equilibrated OMM covered in oil
- Mature cumulus oocyte complexes for 18-24 h (place maturation plate in the back of the incubator)

Prepare media for fertilization

- H-SOF
- 90% Percoll
- SOF-FERT

DAY 0

Preparation of media for fertilization (~2.5 h prior to fertilization)

- ISOLATE (tighten cap and place in warm oven)
 - 1.5 ml 50% over 1.5 ml 90%
- H-SOF (tighten cap and place in warm oven)
 - 1 centrifuge tube per person with 15 ml H-SOF and 1 tube with 10 ml H-SOF (labeled wash)
- SOF-FERT (leave cap loose and place at 38.5°C in 5% CO₂)
 - 35 mm dishes containing 1700 μ l SOF-FERT
 - 1-centrifuge tube with 3 ml SOF-FERT
- PHE (place in oven)
- Warm-up centrifuge canisters (place in oven)
- Plug in citothaw

Matured oocytes: setup for washing and fertilization

- X-Plate with H-SOF
- Searching instrument
- Dissecting microscope
- Heater
- Scissors
- Semen straw plunger
- Inverted microscope
- Small petri-dish
- Rack for tubes (place in front of the heater)
- Slide warmer
- Plastic sterile Pasteur pipets

- Pipette (25 ml)
- Pipet tips

Matured oocytes: washing and fertilization

- Add 5 ml H-SOF to each well of X-plate.
- Transfer about 200 COCs to each well.
- Transfer COCs to 35 mm dish with 1700 μ l/plate SOF-FERT.

Sperm preparation (working in front of the heater)

- Place 1-3 straws of semen in citothaw.
- Layer semen on top of ISOLATE gradient.
- Place ISOLATE tube in a warmed centrifuge canister.
- Centrifuge for 10 min at 1000 x g.
- Collect semen pellet with a Pasteur pipet.
- Place pellet into the 10 ml H-SOF tube.
- Place H-SOF tube into a warmed centrifuge canister.
- Centrifuge for 5 min at 200 x g.
- Pipet off supernatant down to the pellet.
- Add SOF-FERT and determine concentration.

Fertilization

- Add 120 μ l of semen to each fertilization plate containing COCs.
- Add 80 μ l of PHE to each plate.
- Place plates back in the incubator and allow fertilization to proceed for 8-10 h.

Culture media

- Prepare plates with 50 μ l microdrops of SOF-BE1 and cover with mineral oil
- Place plates in the incubator to equilibrate for at least 2 h

Setup for removal of oocytes/embryos from fertilization drop

- Vortexer
- slide warmer
- Hyaluronidase (optional) -warm up
- Sterile dolphin-nose microcentrifuge tubes (and holder)
- Heater (in front of microscope)
- X-Plate
- H-SOF
- Dissecting microscope
- Searching instrument
- Timer

Remove oocytes/embryos from fertilization drop

- Rinse microcentrifuge tube with H-SOF and leave ~30-50 μ l for collection of oocytes/embryos.
- Transfer oocytes/embryos from the fertilization drop to the microcentrifuge tube.
- Vortex the tube for 5 min in front of the heater.
- Use Pasteur pipet to move contents of the microcentrifuge tube to a well of the X-plate.

- Search for cumulus-free oocytes.
- Wash 2X in H-SOF.
- Transfer in groups of up to 30 to the SOF-BE1 microdrops.
- Place culture plate in the back of the incubator.

DAY 3

- Pre-warm stage of inverted microscope and room
- Determine cleavage rate (be quick)

DAY 7-9

- Pre-warm stage of inverted microscope and room
- Collect data on blastocyst development

this page was last updated July 22, 2013
all original material © Rocio Rivera, Peter J. Hansen et al., 2000-2013